

# Mepitel® One



The wound contact layer with Safetac® technology (single sided)


**Safetac**  
TECHNOLOGY

#### Safetac® wound contact layer

- Reduces pain and trauma during dressing changes<sup>1</sup>
- Does not adhere to the moist wound bed, only to dry skin
- Seals the wound margins and reduces risk of maceration<sup>4</sup>


#### Polyurethane net

- Open mesh structure allows for transfer of exudate and application of topical treatment
- High transparency




- Minimizes pain and trauma at dressing changes<sup>1</sup>
- Can remain in place for up to 14 days which allows cost-effective<sup>2,3</sup> and undisturbed wound healing

## How to use Mepitel® One



Clean the wound area. Ensure the peri-wound skin is dry.

Remove the top release film.



Apply Mepitel One to the wound, allowing for at least 2 cm overlap of wounds edges. Larger wounds may require more overlap. Remove second release film during application



Apply outer absorbent dressing and a retention dressing if needed (such as Mextra® Superabsorbent Tubifast®, Tubigrip®, etc.)

## How Mepitel One works

Mepitel One may be left in place for up to 14 days, depending on the condition of the wound, which reduces the necessity for frequent primary dressing changes. The porous structure of Mepitel One allows exudate to pass into an outer absorbent dressing. The Safetac® layer prevents the outer dressing from sticking to the wound and ensures atraumatic dressing changes. The Safetac layer also seals around the wound edges, preventing the exudate from leaking onto the surrounding skin, thus minimizing the risk of maceration.

## Benefits of Mepitel One

- Minimizes pain and trauma at dressing changes<sup>1</sup>
- Safetac layer on wound contact side only (one sided) for easy handling and application
- May remain in place for up to 14 days which allows cost effective<sup>2,3</sup> and undisturbed wound healing
- Enables less frequent dressing changes
- Minimizes the risk of maceration<sup>4</sup>
- Leaves no residue and does not dry out
- Transparent for easy wound inspection during application and wear
- Conforms well to body contours, promoting patient comfort during wear
- Stays in place



Burns



Skin grafts



Skin tears

## Areas of use

- Skin tears
- Surgical incisions
- Traumatic wounds
- Lacerations
- Radiation skin
- Partial and full thickness grafts
- Painful wounds
- Skin abrasions
- Partial thickness burns
- Blistering
- Leg and foot ulcers

It can also be used as protective layer on non-exuding wounds and on areas with fragile skin.

## Wear time

Mepitel One may be left in place for up to 14 days depending on the condition of the wound and surrounding skin, or as indicated by accepted clinical practice (exudate should pass freely through the dressing and the holes should not be blocked).

If saturated, the secondary absorbent dressing should be changed with the Mepitel One left in place.

## Mepitel One ordering information†

Product Code	Size	Pcs/box	Pcs/case	HCPCS
289100	2" x 3" (5 x 7.5 cm)	10	70	A6206
289300	3" x 4" (7.5 x 10 cm)	10	40	A6206
289500	4" x 7" (10 x 18 cm)	10	70	A6207
289700	6.8" x 10" (17 x 25 cm)	5	40	A6208
289750	10.8" x 20" (27.5 x 50 cm)	2	14	A6208
289800	3.75" x 59" (9.5 x 150 cm)	1	9	-

† Packaged sterile in single packs

**References:** 1. White R. et al. Evidence for atraumatic soft silicone wound dressing use. Wounds UK, 2005. 2. Gotschall CS, et al. Prospective, randomized study of the efficacy of Mepitel on children with partial-thickness scalds. Journal of Burn Care & Rehabilitation 1998;19(4):279-283. 3. Rippon M, Davies P, White R, Bosanquet N. Cost implications of using an atraumatic dressing in the treatment of acute wounds. Jo WC, vol 17, No 5 (2008), p. 224-7. 4. Wiberg A.B. et al. Preventing maceration with a soft silicone dressing: in-vitro evaluations. Poster presented at the 3rd Congress of the WUWHS, Toronto, Canada, 2008.

5550 Peachtree Parkway, Suite 500, Norcross, GA 30092 | 1-800-843-8497 | www.molnlycke.us

The Mölnlycke, Mepitel, Tubifast, Tubigrip and Safetac trademarks, names and logo types are registered globally to one or more of the Mölnlycke Health Care Group of Companies. Distributed by Mölnlycke Health Care US, LLC, Norcross, Georgia 30092. © 2018 Mölnlycke Health Care AB. All rights reserved. 1-800-843-8497 MHC-2018-36552

270309 02.18

